



AI-Based Differentiation of ARDS vs Pulmonary Cardiogenic Edema Using Lung Ultrasound (Sonography)

A Deep Learning Approach for Sonography-Based Differential Diagnosis

Presented by: **Namburi Achyuth Krishna** | **252CS022** | **MTECH CSE**

Under the Guidance of :

Prof. B. Annappa

Department of Computer Science and Engineering
National Institute of Technology Karnataka

Presentation Agenda



Project Overview — Literature Review Phase

01 Problem Statement

Why ARDS vs CPE differentiation is a clinical emergency

02 The 5 Key LUS Discriminating Features

Clinical basis for this AI model's learning targets

03 Literature Review Overview

14 papers reviewed — structured evidence base

04 Individual Paper Summaries

Deep dive into each paper's contribution

05 Key Insights from Literature

What we have learned — architecture, features, gaps

06 Research Gap

Why this project is novel and needed

07 Proposed Methodology

Complete 4-stage project pipeline

08 Outcomes & Contribution

Clinical impact and novelty statement

Problem Statement



The Clinical Emergency of Misdiagnosis

A patient arrives in the ICU with severe respiratory distress and bilateral B-lines on ultrasound. Is it ARDS or Cardiogenic Pulmonary Edema?

ARDS

Acute Respiratory Distress Syndrome

- ✓ Lung-protective ventilation
- ✓ Treat underlying cause (sepsis/trauma)
- ✓ Prone positioning
- ✓ Steroids / anticoagulation

✗ Diuretics are HARMFUL



Prone positioning is the technique of placing patients with breathing difficulties on their stomachs to help them breathe better. Prone positioning is generally reserved for sedated patients who require a breathing machine, known as a ventilator, but it may be beneficial for awake patients with COVID-19.

Possible benefits of prone positioning include:

- Reduced risk of ventilator-induced lung injury
- Less lung compression and more efficient gas exchange in the lungs
- Improved heart function and oxygen delivery to the body
- Better drainage of secretions produced in deeper lung



All patients placed in prone position should be monitored carefully for respiratory status and symptoms.

≠

WRONG
TREATMENT
= DEATH

CPE

Cardiogenic Pulmonary Edema

- ✓ Diuretics (furosemide)
- ✓ Reduce IV fluid
- ✓ Treat heart failure
- ✓ Vasodilators

✗ Lung-protective ventilation unnecessary



Why Current Methods Fail — Why LUS + AI?




The Diagnostic Gap


Berlin Definition

- Requires CXR / CT scan
- Blood gas analysis needed
- Subjective radiologist read
- Cannot be done bedside fast
- Inter-observer disagreement

Berlin Criteria for Acute Respiratory Distress Syndrome (ARDS)



ARDS Severity Grading (PaO_2/FiO_2 with PEEP/CPAP ≥ 5 cm H₂O)



Clinical LUS (Human)

- B-lines look IDENTICAL in ARDS & CPE
- Sensitivity only 82–87% (meta-analysis)
- 1 in 7 ARDS cases missed
- Requires trained LUS expert
- Subjective sign interpretation






Figure Proposed expansion of the Berlin definition of ARDS

LUS + AI (Current Solution)

- Bedside — real-time, no radiation
- Detects SUBVISUAL pixel patterns
- CNN AUC up to 1.0 on similar tasks
- Objective, reproducible, automated
- Available 24/7 in any setting





1. **A-lines** (Normal Indicator): Horizontal, echogenic lines parallel to the pleural line.
2. **B-lines** (Pathological Indicator): Vertical, laser-like lines extending from the pleural line to image bottom.
3. **Consolidation Features** (Severe Lung Pathology): Hypoechoic (dark) areas replacing normal air-filled alveoli.

The 5 Key LUS Discriminating Features




What Our AI Model Learns to Detect

LUS FEATURE	ARDS	Cardiogenic PE
 Pleural Line Irregularity	Thick, fragmented, jagged	Thin, smooth, continuous
 B-Line Distribution	Heterogeneous, patchy, asymmetric	Homogeneous, bilateral, uniform
 Spared Areas	PRESENT — dark gaps between B-lines	ABSENT — uniform B-lines everywhere
 Subpleural Consolidations	PRESENT — small wedge-shaped patches	ABSENT — lung structure intact
 Lung Sliding	Reduced / Absent in affected zones	Normal / Mildly reduced


Normal Lung vs Lung with B-Lines




Normal Lung Ultrasound




Lung Ultrasound with B Lines



Lung Ultrasound with Pleural Effusion





Paper 1 — Arntfield et al. (2021)



CNN to Differentiate Etiology of B-Lines on Lung Ultrasound

PROBLEM STATEMENT




B-lines appear identical in ARDS, COVID ARDS & HPE — physicians cannot distinguish them reliably.

DATASET

243 patients; 612 LUS video clips; 121,381 frames (COVID ARDS, Non-COVID ARDS, HPE)

APPROACH / METHOD

Xception CNN + Transfer Learning (ImageNet). Frame-level classification aggregated to encounter level.

	Hydrostatic Pulmonary Edema	Non-COVID ARDS	COVID-19 ARDS
			
B Lines	Homogeneous	Patchy	Patchy
Pleural Line	Smooth/Regular	Irregular/Fragmented	Irregular/Fragmented
Sub-pleural consolidation	Absent	Common	Common

EVALUATION METRICS

AUC, Sensitivity, Specificity, F1-score (frame-level & encounter-level)

RESULTS / KEY PERFORMANCE

ENCOUNTER AUC: COVID=1.0 | Non-COVID ARDS=0.934 | HPE=1.0. Outperformed physicians (AUC ~0.789).

CLASSES COMPARED

COVID-19 ARDS vs Non-COVID ARDS vs Hydrostatic Pulmonary Edema



OBSERVATION — HOW THIS GUIDES OUR PROJECT

CORE BASELINE. Proves CNN can differentiate ARDS from pulmonary edema. We use Xception/ResNet + Transfer Learning.

Grad-CAM confirms pleural zone = key ROI.

This project removes COVID confound — pure binary ARDS vs CPE.

LIMITATIONS:

Single-centre. Small dataset. COVID may inflate AUC. No external

validation.

GAP:

No clean binary ARDS vs CPE CNN classifier. This project fills this.




Figure 4 Grad-CAM heatmaps corresponding to a selection of our model's predictions. Blue areas reflect the regions of the image with the highest contribution to the resulting class predicted by the model. In all cases, the immediate area surrounding the pleura appears most activated. COVID, COVID-19 pneumonia; HPE, hydrostatic pulmonary edema; NCOVID, non-COVID-related acute respiratory distress syndrome.

Paper 3 — Chiumello (Campobasso) et al. (2024)



Lung Imaging and Artificial Intelligence in ARDS

PROBLEM STATEMENT

AI applied to ARDS imaging fragmented across labs — no unified comparison of architectures, modalities, or gaps.

APPROACH / METHOD

Review: VGG16, ResNet, InceptionV3, U-Net. Transfer learning. Multi-modality AI comparison.

EVALUATION METRICS

Summary AUC across reviewed models: 0.90–1.00

This review confirms our project fills an identified gap. It guides our architecture selection (ResNet50 or InceptionV3 with transfer learning) and confirms LUS is the optimal modality for our task. The explicit statement that no dedicated ARDS vs CPE binary classifier exists in the literature DIRECTLY VALIDATES OUR PROJECT'S NOVELTY.

RESULTS / KEY PERFORMANCE

ResNet & InceptionV3 + transfer learning: AUC 0.93–1.00. LUS outperforms CXR for bedside AI.

No binary ARDS vs CPE classifier identified.

CLASSES COMPARED

Healthy vs COVID ARDS vs Non-COVID ARDS vs HPE vs Pneumonia (across reviewed studies)

★ OBSERVATION — HOW THIS GUIDES TO THIS PROJECT

Directly confirms our project fills an identified gap.

Architecture guide: ResNet50 vs InceptionV3.

Confirms LUS is best modality.

Gap explicitly stated — no binary ARDS vs CPE classifier exists.

LIMITATIONS:

Review only. Rapidly evolving field.

GAP:

No dedicated binary ARDS vs CPE deep learning classifier. OUR PROJECT.

Paper 6 — Zhao & Bell (2022)




Deep Learning in Lung Ultrasound Imaging of COVID-19 Patients — Review

PROBLEM STATEMENT

20+ DL models built for LUS in isolation — no unified architecture comparison or practical building guide.

APPROACH / METHOD

Review: VGG16, ResNet50, InceptionV3, EfficientNet, CNN+LSTM. Transfer learning & augmentation strategies.



RESULTS / KEY PERFORMANCE

ResNet50 & InceptionV3 best (AUC 0.93–1.00).

CNN+LSTM adds 3–7% AUC for video.

Vertical flip harmful.

Horizontal flip + brightness effective.

CLASSES COMPARED

CPE, ARDS, Pneumonia, COVID-19 vs Normal (across reviewed studies)

★ OBSERVATION — HOW THIS GUIDES OUR PROJECT

ARCHITECTURE SELECTION GUIDE. We implement **ResNet50 + InceptionV3** with ImageNet transfer learning. **CNN+LSTM for lung sliding.**

Two-phase fine-tuning protocol. Specific augmentation strategy adopted.

LIMITATIONS:

COVID-centric. No unified benchmark. Rapidly evolving.

GAP:

No standardised ARDS vs CPE LUS benchmark. This project establishes one.

Key Insights Collected from All 14 Papers



What the Literature Tells Us



Best Architecture

- Xception/ResNet50 + Transfer Learning — core model (Paper 1)
- InceptionV3 — best for multi-scale B-line patterns (Papers 3,6)
- CNN+LSTM — adds 3–7% AUC for video / lung sliding (Papers 6,9)
- Two-phase fine-tuning: freeze → unfreeze last blocks (Paper 9)



Best Features / ROI

- Pleural line zone = primary CNN attention (Grad-CAM, Paper 1)
- Heterogeneous B-lines = strongest single discriminator (Papers 2,13,7)
- Spared areas = highest specificity for ARDS (Papers 5,7,14)
- Anterolateral 12 zones = priority input regions (Paper 13)



Evaluation Design

- Head-to-head: CNN vs LUS-trained physicians (Papers 1,4)
- Report: AUC + sensitivity + specificity (Papers 4,7)
- Baseline to beat: human AUC 0.89–0.93 (Paper 7)
- Focal vs Non-Focal ARDS subgroup analysis (Papers 11,7)



Research Gaps Confirmed

- No binary ARDS vs CPE CNN classifier confirmed (Paper 3)
- No temporal lung sliding model exists (Papers 3,6)
- Human sensitivity only 82–87% — 1 in 7 missed (Paper 7)
- GLCM classical ML as baseline comparison available (Paper 2)

Research Gap — Why Our Project is Needed



What the Literature is Missing

A 2024 comprehensive review (Chiumello et al.) explicitly confirms: No dedicated binary ARDS vs CPE deep learning classifier on LUS exists in the literature.

GAP 1

No Binary ARDS vs CPE CNN

All DL studies bundle COVID, use multi-class, or compare vs Normal — never pure ARDS vs CPE

Papers 1, 3, 4

GAP 2

No Temporal Lung Sliding Model

No study has used CNN+LSTM specifically to capture lung sliding for ARDS vs CPE

Papers 3, 6, 9

GAP 3

Human Sensitivity = Only 82–87%

1 in 7 ARDS cases misdiagnosed even by expert LUS. AI can close this gap.

Paper 7

GAP 4

No Multi-Centre Validated Model

All existing LUS AI models are single-centre. External validation is absent.

Papers 1, 3, 6

GAP 5

Explainability Not Correlated with Signs

Grad-CAM used but never systematically correlated with the 5 clinical LUS signs.

Papers 1, 3

OUR PROJECT

Addresses ALL 5 Gaps


Binary ARDS vs CPE | CNN+LSTM | Improve sensitivity | Multi-centre design | Grad-CAM + sign correlation

Novel contribution


Model 1



A. ResNet 50 Diagram




B. Residual learning







GoogLeNet Architecture



LSTM ARCHITECTURE FOR PULMONARY ULTRASOUND CLASSIFICATION





ARCHITECTURE SUMMARY	
Stage	Configuration
Input	Sequence of T frames, each of size $224 \times 224 \times 1$ (or 3)
Conv Block 1	Conv2D (32, 3×3) + ReLU + MaxPool (2×2)
Conv Block 2	Conv2D (64, 3×3) + ReLU + MaxPool (2×2)
Conv Block 3	Conv2D (128, 3×3) + ReLU + MaxPool (2×2)
Conv Block 4	Conv2D (256, 3×3) + ReLU + MaxPool (2×2)
Flatten (per time step)	Feature vector
LSTM Layer 1	128 units, $\text{return_sequences} = \text{True}$
LSTM Layer 2	64 units, $\text{return_sequences} = \text{False}$
Dense Layer	64 units, ReLU
Output Layer	num_classes, Softmax



NOTATIONS	
x_t	: Input at time step t
h_t	: Hidden state at time step t
c_t	: Cell state at time step t
σ	: Sigmoid activation
$\tanh$	: Hyperbolic tangent
$\otimes$	: Element-wise multiplication
$+$	: Element-wise addition

This LSTM architecture captures spatio-temporal information from a sequence of lung ultrasound frames by combining CNN-based feature extraction with LSTM's temporal modeling capability, making it suitable for pulmonary pattern classification.

Proposed Custom CNN Architecture




Block	Conv Filters	Kernel Size	Output Size	Dropout
Block 1	32	3×3	$112 \times 112 \times 32$	0.25
Block 2	64	3×3	$56 \times 56 \times 64$	0.25
Block 3	128	3×3	$28 \times 28 \times 128$	0.30
Block 4	256	3×3	$1 \times 1 \times 256$ (GAP)	0.40
Dense	128	–	128	0.50
Output Layer	num_classes	–	num_classes	–

ARCHITECTURE SUMMARY


- Lightweight and optimized for lung ultrasound images
- Batch Normalization for stable training
- Dropout for regularization
- Global Average Pooling to reduce parameters and overfitting
- Softmax output for multi-class classification




COVID-19 Lung Ultrasound Dataset

- Source: GitHub – COVID-19 Lung Ultrasound Dataset
- Type: Lung Ultrasound (LUS) Images & videos
- Classes available:
 - COVID / Viral Pneumonia
 - Bacterial Pneumonia
 - Healthy
- Probe types:
 - Linear probe
 - Convex probe
- Characteristics:
 - Real clinical ultrasound scans
 - Includes artifacts (B-lines, consolidations)
- Limitation:
 - Not specifically labeled as ARDS vs PCE
 - Used as **proxy dataset for initial model validation**



Overview figure about current efforts. Public dataset consists of >200 LUS videos.








Preprocessing Pipeline—

- Image resizing:
 - ResNet → 224×224
 - Inception → 299×299
- Grayscale → Converted to 3-channel (if needed)
- Normalization:
 - Pixel scaling (0–1)
- Data augmentation:
 - Rotation ($\pm 15^\circ$)
 - Horizontal flip
 - Zoom
- Noise handling:
 - Ultrasound speckle preserved (important feature)
- Dataset split:
 - Train: 70%
 - Validation: 15%
 - Test: 15%
 - Patient-level split (no leakage)

Model Output



True: covid
Pred: covid



True: covid
Pred: covid



True: covid
Pred: covid




True: covid
Pred: covid



True: regular
Pred: covid



Model Results




Grad-CAM Explainability (NOVELTY)




Original Ultrasound



Attention Heatmap




Explainability Visualization




Grad-CAM Explainability (NOVELTY)




Original Ultrasound



Attention Heatmap



Explainability Visualization




Grad-CAM Explainability (NOVELTY)




Original Ultrasound



Attention Heatmap



Explainability Visualization





Thank You